

# COPPER

## (Colorimetric Test with Dibrom-PAESA)

Cat.No	Package Size
162 000	2 x 50 mL + Standard
162 016	7 x 10 mL + Standard

### PRINCIPLE:

Copper forms with 4-(3,5-dibromo-2-pyridylazo)-N-ethyl-N-sulfopropylaniline a chelate complex. The increase of absorbance of this complex can be measured and is proportional to the concentration of total copper in the sample.

### REAGENTS (ready for use)

Monoreagent		
acetate buffer pH 5.0		0.2
mol/l		
4-(3,5-dibromo-2-pyridylazo)- N-ethyl-N-sulfopropylaniline		0.02
mmol/l		
<b>Standard</b>	100 µg/dl	= 15,7
µmol/l		

### STORAGE AND STABILITY:

The sealed reagent is stable up to the indicated expiry date if stored at 2 - 25°C.

#### NOTE:

*If stored at 2°- 8°C precipitation may occur. In this case store the reagent at for about 2 hours over 20°C and mix until the reagent is clear.*

### SAMPLE MATERIAL:

Serum, Plasma

### QUALITY CONTROL:

Use Greiner's control sera with **Copper values**, determined by this method or by atomic absorbance.

### REFERENCE VALUES:

70 - 153 µg/dl in serum (11 -24 µmol/l)

### LINEARITY:

Up to 500 µg/dl (78.65 µmol/l)

### INTERFERENCES :

By in-vitro tests it was shown that only Nickel (Ni) does interfere heavily, but Nickel is hardly present in normal human samples.

### ANALYTICAL PROCEDURE

Wavelength : 580 nm  
Light path: 1 cm  
Temperature : 37°C

Pipette into cuvettes:

	Standard	Sample
<b>Monoreagent</b>	1000 µl	1000 µl
<b>Serum or plasma</b>	-	50 µl
<b>Standard</b>	50 µl	-

Mix and incubate for 5 minutes at 37°C. Measure the absorbance of the sample  $A_S$  and of the standard  $A_{STD}$  against the reagent blank

$$A_{RBL} \\ \Delta A_S = A_S - A_{RBL} \\ \Delta A_{STD} = A_{STD} - A_{RBL}$$

### CALCULATION:

$$\mu\text{g/dl copper} = \frac{\Delta A_S}{\Delta A_{STD}} \times 100$$

$$\mu\text{mol/l copper} = \frac{\Delta A_S}{\Delta A_{STD}} \times 15,7$$

### LITERATURE:

Abe A., Yamashita S., Noma A., Clin. Chem., 552-554- 35 (1989)

### SYMBOLS USED

 For *in vitro* diagnostic medical use

 Batch Code

 Use by

 Temperature limitation