


# OXALATE

## (Quantitative Enzymatic Color Test)

Cat.No	Package Size		
	R1	R2	Sample Diluent
175 080	10 x 10 ml	5 x 2 ml + 2 x 10 ml Standard	5 x 100 ml
175 081	2 x 10 ml	1 x 2 ml + 4 ml Standard	1 x 100 ml
175 082	Sample Purifier Tubes 20 pieces		

### METHOD / TESTPRINCIPLE

Enzymatic colorimetric Test :



### REAGENTS COMPOSITION (concentration in test):

<b>R1 :</b>	DMAB	3.2 mmol/L
	MBTH	0.22 mmol/L
	Buffer	pH 3.1
<b>R2 :</b>	Oxalate Oxidase	3.0 kU/L
	Peroxidase	100 kU/L
<b>Sample Diluent :</b>		
	EDTA	10 mmol/L
	Buffer	pH 7.6
<b>Standard :</b>		
	Oxalate	0.5 mmol/L

### Additional Reagent (Cat. No 175 082)

not included in the kit :

**Sample Purifier tubes** Activated Charcoal

### STABILITY OF REAGENTS

When stored at 2-8° C and protected from light and contamination, the reagents are stable up to the expiry date printed on the labels. Do not freeze !

### PREPARATION AND STABILITY OF WORKING REAGENTS

**R1 :** Reconstitute Substrate with 10 ml of dist. water.

**R2 :** Reconstitute Start Reagent with 2 ml of dist. water .

Stability of R1 : 2 weeks at 2-8° C

R2 : 5 weeks at 2-8° C

**Sample Diluent** is ready to use

### SAMPLES

Urine

### REFERENCE VALUES

	mg/24 hrs	mmol/24 hrs
Adult Males	7 - 44	0.08 - 0.49
Adult Females	4 - 31	0.04 - 0.32
Children	13 - 38	0.14 - 0.42

### PROCEDURE

#### Sample preparation:

1. Set up a series of labelled tubes for urine samples and controls.
2. Pipette 5 ml or any suitable volume of urine samples and controls into **Sample Purifier Tubes**.
3. Add equal volume (as in step 2) of *sample diluent* into each tube and mix for approximately 5 min by intermittent mixing (rotator mixer is recommended).
4. Check the pH, it should be between 5.0 and 7.0 , otherwise adjust the pH using 1N hydrochloric acid and/or 1N sodium hydroxide
5. Centrifuge tubes for 5 min at 2000 rpm (1500 g) or use filtration.  
Determine oxalate concentration in the supernatants or filtrates as described below.

#### Determination of Oxalate

1. Warm oxalate reagents to assay temperature (preferably 25° C or 30° C or 37° C).
2. Label tubes for *reagent blank*, *standard*, *controls* and *urine samples*
3. Pipette 1 ml **R1** into each tube.
4. Pipette 50 µl of *supernatants* (or *filtrates*) , to respective tubes. Add 50 µl *deionized water* to reagent blank tube and 50 µl *standard* to tube labelled standard.
5. Pipette 0.1 ml of **R2** into each tube and immediately mix by gentle inversion.
6. Incubate the tubes at desired temperature for 5 minutes.
7. Read absorbances (**A**) of *blank*, *standard*, *controls* and *urine samples* at 590 nm.
8. Determine the corrected absorbances (**ΔA**) of standard, controls and sample by subtracting reagent blank absorbance from the *absorbance* readings of *standard*, *controls* and *urine samples*.

#### Calculation

$$\text{Oxalate (mmol/L)} = \frac{\Delta A \text{ Sample}}{\Delta A \text{ Standard}} \times 0.5 \times 2$$

0.5 = Concentration (mmol/L) of oxalate in standard

2 = Dilution factor

#### Quantity of Oxalate excreted during 24 hrs =

$$\text{Oxalate (mmol/L)} \times \text{Volume of urine (Litres)}$$

### QUALITY CONTROL

For quality control use adequate control materials.