

PHI

Kinetic UV-Test for the determination of Phosphohexose Isomerase in serum and plasma

Cat. No	Package Size
178 080	R1 (Buffer) 4 x 6 ml / R2 (Lyophilisate) for 4 x 6 ml

Method / Testprinciple

Kinetic UV-Test according to Bueding and Mackinnon:

Fructose-6-phosphate $\xrightarrow{\text{PHI}}$ Glucose-6-phosphate

Glucose-6-phosphate + NAD⁺
 $\xrightarrow{\text{G6P-DH}}$ 6-Phosphogluconate + NADH₂⁺

Reagents

Constituents (concentrations in the test)

R1: TRIS-buffer (pH 8.5)	100 mmol/l
R2: Fructose-6-phosphate	3 mmol/l
Nicotinamide-adenin-dinucleotide (NAD)	3 mmol/l
Glucose-6-phosphate-dehydrogenase (G6P-DH)	1 U/ml

Storage and stability

The reagents are stable up to the expiry date printed on the labels, when stored at 2-8° C and protected from light.

After reconstitution the reagent solution is stable at least for 2 weeks at 2 – 8°C, when contamination is avoided.

Precautions

- For in vitro diagnostic use only.
- The reagent contains < 0,95 g/L sodium azide.
- Avoid contact with skin and/or mucous membranes.

Waste

Treat according to the local legal regulations

Preparation

Reagent R :

Reconstitute 1 bottle of lyophilisate R2 („for 6 ml“) with 1 bottle of buffer R1.

R is after 15 min ready for use.

Stability in tightly closed bottle and without contamination: :

- 4 days at 2 – 8 °C
- 1 day at 18 – 25 °C

Samples

Serum, plasma, liquor

Stability of PHI in serum 2 weeks at 2 - 8 °C

*3 months at < - 25 °C
(Freeze only once!)*

Procedure

Wavelength	340 nm, Hg 334 nm, Hg 365 nm
Cuvette	1 cm
Temperature	37°C (25 °C)
Measure	against air (increasing absorbance)

Pipette into cuvettes

	Semi-Micro	Micro
Reagent R	1000 µl	500 µl
Sample	50 µl	20 µl

Mix well, incubate 1 min at constant temperature, read absorbance A and immediately start stop watch. After exactly 1, 2 and 3 min read again. Calculate ΔA/min

Calculation

$$\Delta A/\text{min} \times \text{Factor} = \text{PHI-Activity [U/I]}$$

Factors for 37 °C (25 °C)			
Wavelength	Semi-Micro		Micro
Hg 334 nm	6150	(3398)	7615 (4207)
340 nm	6033	(3333)	7470 (4127)
Hg 365 nm	1117	(6176)	13841 (7647)

Applications for automated systems are available on request

Calibration and Quality Control

For the calibration of automated analyzers Greiner multicalibrator Unical-M is recommended, for quality control use Greiner normal and abnormal controls Unitrol-I and Unitrol-II respectively the PHI special control