

# PYRUVATE

## (Quantitative Enzymatic UV-Test)

Cat.No.	Package Size		
	R1a	R1b	R2
180 000	2 x 50 mL	1 x 5.0 mL	1 x 5.0 mL
	Standard 1 x 20 mL		

### METHOD / TESTPRINCIPLE

Enzymatic. UV-Test .

Endpoint determination of pyruvate in blood :

In the presence of an excess of NADH pyruvate is converted to lactate.

The reduction of the absorbance =  $\Delta A$  ,at 340 nm, due to the oxidation of NADH to NAD<sup>+</sup> , is a measure of the amount of pyruvate originally present :



*LDH = Lactatedehydrogenase*

### REAGENTS COMPOSITION

R1a :	(Tris buffer, pH 7.20)	1.50 mol/L
R1b :	NADH	10.0 mmol/L
Start-Reagent R2 :	LDH	1.50 kU/mL
Standardsolution	Pyruvate	4.00 mg/dL

**Additional Reagent** ( not provided with the kit )

**Cat.No.502 001 (2 x 500ml) Deproteinization Reagent (Perchloric Acid 0.6m)**

### PRECAUTIONS

- For in vitro diagnostic use only.
- Reagents contain < 0,95 g/ L sodium azide. Avoid contact with skin and/or mucous membranes
- Avoid contamination by using clean laboratory material (pipettes, plastic vials for analyzers, ...).
- Discard cloudy i.e. deteriorated reagent.

### STABILITY OF REAGENTS

When stored at 2-8° C and protected from light, the reagents are stable until the expiry date stated on the labels.

### PREPARATION AND STABILITY OF WORKING REAGENTS

Reagents (R1a ,R1b and reagent R2) are ready for use

### SAMPLE PREPARATION

Pipet 2,0 mL of **freshly drawn blood** into a centrifugation tube containing 4 mL of cold **0.6 m perchloric acid**. Vortex for about 30 seconds. Keep the blood precipitate mixture for about 5 min in the cold to assure complete protein precipitation. Centrifuge 10 min at approximately 1500 x g. The protein free supernatant is ready for use.

The **Standardsolution** has to be diluted with perchloric acid, in the same ratio as the sample.

### MANUAL PROCEDURE

Wavelength : 340 nm (334-365)

Temperature : 30°C, 37°C

Cuvette : 1 cm light path  
(Measure against water)

	Supernatant Sample	Diluted Standard
Buffer Reagent R1a	1.00 mL	1.00 mL
Supernatant Sample	2.00 mL	-
Diluted Standard	-	2.00 mL
Mix and add		
Reagent R1b	50 µL	50 µL
Mix and incubate for approx. 5 min . Pour into cuvette, measure initial absorbance A <sub>1</sub> .		
Start-Reagent R2	50 µL	50 µL
Mix, incubate for approx. 5 min and measure absorbance A <sub>2</sub>		

$$\Delta A = A_2 - A_1 \text{ Sample and Standard}$$

### CALCULATION

- with Factor :

$$\text{Pyruvate (mg/dL)} = \Delta A_{340} \times 6,37$$

- with Standard :

$$\text{Pyruvate (mg/dL)} = \Delta A_{\text{sample}} \times 4.00 / \Delta A_{\text{Standard}}$$

### REFERENCE VALUES (37° C) :

**0,3 – 0,7 mg/dL**

*Note: It is recommended for each laboratory to establish and maintain its own reference values. The given data are only an indication.*

### QUALITY CONTROL

For quality control use adequate control materials, available from Greiner .